

creates 3 targets in a linkmodelbased rotation coordinate system.

target 1 = origin

target 2 = origin - (0,0,marker_distance)

target 3 = origin + (0,marker_distance,0)

```
create_targets_for_linkmodelbased_rotations
/signal_folders=
! /marker_distance=0.05
;
```

[back to pipeline commands reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:create_targets_for_linkmodelbased_rotations&rev=1718804219

Last update: 2024/06/19 13:36

