

Enter Password to View Folder Data

The `Enter_Password_To_View_Folder_Data` allows a user to decrypt a password-protected data tree folder to view its contents.

Pipeline Command

```
Enter_Password_To_View_Folder_Data
/SIGNAL_TYPES=
/SIGNAL_FOLDER=ORIGINAL
/PASSWORD=
;
```

Command Parameters

Parameter	Description
Signal Types	The type of signal whose folder will be accessed.
Signal Folder	The signal folder to access.
Password	The folder's password.

See Also

Related pipeline commands include:

- [Add_Password_Protection_To_Folder](#)
- [Remove_Password_Protection_From_Folder](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:enter_password_to_view_folder_data

Last update: 2026/09/10 19:12

