

fixes a mac eva7 c3d file. if the file was modified in the c3d editor use original\_eva7=false.

```
fix_mac_eva7_c3dfile
! /file_name=
! /suffix=_
! /original_eva7=true
;
```

[back to pipeline commands reference](#)

From:  
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file\\_commands:fix\\_mac\\_eva7\\_c3dfile&rev=1718804227](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:fix_mac_eva7_c3dfile&rev=1718804227)

Last update: **2024/06/19 13:37**

