

imports analog signals from an anb formatted file into a c3d file.

```
import_analog_anb_signals  
! /replace_existing_signals=true  
! /delete_all_existing_signals=false  
! /update_c3d_file=false  
;
```

the filename must be the same; only the extension is different.
for example:

file.c3d will import signals from **file.anb**

[back to pipeline commands reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_anb_signals&rev=1718804227

Last update: 2024/06/19 13:37

