

Imports analog data from a Bioware (kistler) v3 file (.dat) into a C3D file.

```
Import_Analog_Bioware_Signals  
! /FILE_NAME=  
;
```

[Back to Pipeline Commands Reference](#)

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_bioware_signals&rev=1719432323

Last update: **2024/06/26 20:05**

