

# Import\_Analog\_Bioware\_Signals

Imports analog data from a Bioware (kistler) v3 file (.dat) into a C3D file.

```
Import_Analog_Bioware_Signals
! /FILE_NAME=
;
```

[Back to Pipeline Commands Reference](#)

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file\\_commands:import\\_analog\\_bioware\\_signals&rev=1720027709](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_bioware_signals&rev=1720027709)

Last update: **2024/07/03 17:28**

