

Import_Analog_Bioware_Signals

Imports analog data from a Bioware (kistler) v3 file (.dat) into a C3D file.

```
Import_Analog_Bioware_Signals
! /FILE_NAME=
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_bioware_signals&rev=1721158125

Last update: **2024/07/16 19:28**

