

Import Analog Bioware Signals

Imports analog data from a Bioware (kistler) v3 file (.dat) into a C3D file.

```
Import_Analog_Bioware_Signals
! /FILE_NAME=
! /SET_PROMPT=Open data file
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_bioware_signals&rev=1741027056

Last update: 2025/03/03 18:37

