

# Import Analog Signals From DASYLAB Ascii File

Imports analog signals from a DASYLab ASCII formatted file into a C3D file.

```
Import_Analog_Signals_From_DASYLab_Ascii_File  
! /REPLACE_EXISTING_SIGNALS=TRUE  
! /FOLDER=  
;
```

Note that the name of the analog file should be similar to the C3D file.  
For example: **test.c3d** and **test\_a.txt**

[Back to Pipeline Commands Reference](#)

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file\\_commands:import\\_analog\\_signals\\_from\\_dasylib\\_ascii\\_file](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_signals_from_dasylib_ascii_file)

Last update: 2024/07/17 15:45

