

Import_Analog_SON_Signals

Imports analog data from a SON data file (.smr/spike2) into a C3D file.

```
Import_Analog_SON_Signals
! /FILE_NAME=
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_analog_son_signals&rev=1720027714

Last update: **2024/07/03 17:28**

