

imports data from a foot pressure sensor.

```
import_pedar_data  
/right_def_file=  
/right_data_file=  
/c3d_file=  
/left_def_file=  
/left_data_file=  
/shoe_cal_file=  
! /insole_target_numbers=  
;
```

[back to pipeline commands reference](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_pedar_data&rev=1718801430

Last update: **2024/06/19 12:50**

