

Import_Pedar_Data

Imports data from a foot pressure sensor.

```
Import_Pedar_Data  
/RIGHT_DEF_FILE=  
/RIGHT_DATA_FILE=  
/C3D_FILE=  
/LEFT_DEF_FILE=  
/LEFT_DATA_FILE=  
/SHOE_CAL_FILE=  
! /INSOLE_TARGET_NUMBERS=  
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_pedar_data&rev=1720790978

Last update: **2024/07/12 13:29**

