

Imports data from University of Surrey's Force Platform channels into a C3D file.

```
Import_Surrey_Force_Platform_Signals  
! /OVERWRITE_SIGNALS=TRUE  
;
```

[Back to Pipeline Commands Reference](#)

From:  
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file\\_commands:import\\_surrey\\_force\\_platform\\_signals&rev=1718805613](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_surrey_force_platform_signals&rev=1718805613)

Last update: **2024/06/19 14:00**

