

Import_Surrey_Force_Platform_Signals

Imports data from University of Surrey's Force Platform channels into a C3D file.

```
Import_Surrey_Force_Platform_Signals  
! /OVERWRITE_SIGNALS=TRUE  
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:import_surrey_force_platform_signals&rev=1720028443

Last update: **2024/07/03 17:40**

