

this command is not available for common use.

creates linkmodelbased joint rotations based on psu imu data.

target 1 = origin

target 2 = origin - (0,0,marker_distance)

target 3 = origin + (0,marker_distance,0)

```
joint_rotations_from_psu_imu
/motion_file=
/floor_file=
/subject_cal_file=
/c3d_file=
;
```

[back to pipeline commands reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:joint_rotations_from_psu_imu&rev=1718801432

Last update: **2024/06/19 12:50**

