

merges a set of ndi aurora files into a visual3d ascii file.

```
merge_aurora_files
! /input_file_names=
! /output_file_name=
! /frame_rate=1000
! /units=cm
;
```

wildcards (*) can be used for the input_file_names.

[back to pipeline commands reference](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:merge_aurora_files&rev=1718804236

Last update: **2024/06/19 13:37**

