

Modifies the C3D MANUFACTURER Group in a C3D File. Optionally, overwrite the original file on disk.

[Modify_C3D_Manufacturer.png](#)

```
Modify_C3D_Manufacturer
! /COMPANY=
! /SOFTWARE=
! /VERSION_LABEL=
! /VERSION=
! /OVERWRITE_C3D_FILE=FALSE
;
```

[Back to Pipeline Commands Reference](#)

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:modify_c3d_manufacturer&rev=1718648376

Last update: **2024/06/17 18:19**

