

Modify_C3D_Manufacturer

Modifies the C3D MANUFACTURER Group in a C3D File. Optionally, overwrite the original file on disk.



```
Modify_C3D_Manufacturer
! /COMPANY=
! /SOFTWARE=
! /VERSION_LABEL=
! /VERSION=
! /OVERWRITE_C3D_FILE=FALSE
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:modify_c3d_manufacturer&rev=1720793024

Last update: 2024/07/12 14:03

