

Modify_C3D_Manufacturer

Modifies the C3D MANUFACTURER Group in a C3D File. Optionally, overwrite the original file on disk.

Identify the Manufacturer of This C3D File

The C3D file format was established many years ago (see www.c3d.org). Over time there have been changes to the kinds of data stored in the file, and there have been unfortunate mistakes in some versions of software (including Visual3D).

Consequently, it is useful to store software names (and sometimes versions) in the C3D file. In some cases only the Manufacturer name is sufficient information. In some cases the anomaly may be specific to only one release version of one software product.

This dialog will only request sufficient information, so please fill in all information that is relevant.

Manufacturer

Get Current Parameters

OK

Cancel

Modify_C3D_Manufacturer

```
! /COMPANY=  
! /SOFTWARE=  
! /VERSION_LABEL=  
! /VERSION=  
! /OVERWRITE_C3D_FILE=FALSE  
;
```

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:modify_c3d_manufacturer&rev=1721229056

Last update: 2024/07/17 15:10

