

Modify the parameters in the SUBJECTS group in the c3d file.

Modify_C3D_Subjects_Parameters

```
! /USED=  
! /IS_STATIC=  
! /USES_PREFIXES=  
! /NAMES=  
! /LABEL_PREFIXES=  
! /MARKER_SETS=  
! /DISPLAY_SETS=  
! /MODELS=  
! /MODEL_PARAMS=  
! /OVERWRITE_C3D_FILE=FALSE  
;
```

[c3dSubjectDlg.jpg](#)

[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:modify_c3d_subjects_parameters&rev=1718648377

Last update: 2024/06/17 18:19

