

Modify the parameters in the SUBJECTS group in the c3d file.

#### Modify\_C3D\_Subjects\_Parameters

```
! /USED=  
! /IS_STATIC=  
! /USES_PREFIXES=  
! /NAMES=  
! /LABEL_PREFIXES=  
! /MARKER_SETS=  
! /DISPLAY_SETS=  
! /MODELS=  
! /MODEL_PARAMS=  
! /OVERWRITE_C3D_FILE=FALSE  
;
```



[Back to Pipeline Commands Reference](#)

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file\\_commands:modify\\_c3d\\_subjects\\_parameters&rev=1719432334](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:modify_c3d_subjects_parameters&rev=1719432334)

Last update: **2024/06/26 20:05**

