

Modify_C3D_Subjects_Parameters

Modify the parameters in the SUBJECTS group in the c3d file.

```
Modify_C3D_Subjects_Parameters
```

```
! /USED=  
! /IS_STATIC=  
! /USES_PREFIXES=  
! /NAMES=  
! /LABEL_PREFIXES=  
! /MARKER_SETS=  
! /DISPLAY_SETS=  
! /MODELS=  
! /MODEL_PARAMS=  
! /OVERWRITE_C3D_FILE=FALSE  
;
```



[Back to Pipeline Commands Reference](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:file_commands:modify_c3d_subjects_parameters&rev=1721158141

Last update: **2024/07/16 19:29**

