

in version 5 we introduced the parameter **component_sequence**

this replaced the legacy parameter **components**

there were two reasons behind this change.

1. it allows the user to specify the order of the components

/component_sequence= x + z + y 2. it allows the user to specify more than 4 components

/component_sequence= 0 + 2 + 4 + 6

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:general_information:component_sequence&rev=1718804252

Last update: **2024/06/19 13:37**

