

Signal Folder

This is a common parameter in pipeline commands, allowing the user to specify which [data tree folder](#) contains the signals to be evaluated as input. Signal folder is the second (middle) level of the [data tree](#) hierarchy.

The parameter will appear in commands as `/Signal_Folder` or `/SIGNAL_FOLDER`.

For example:

- ORIGINAL
- PROCESSED
- User-named folders, etc.

See Also

The `Signal_Type` parameter is often accompanied by two other parameters:

- [Signal_Types](#)
- [Signal_Names](#)

Together these three parameters represent the application-wide [data tree syntax](#).

From:
<https://www.wiki.has-motion.com/> - **Wiki Documentation Site**

Permanent link:
https://www.wiki.has-motion.com/doku.php?id=visual3d:documentation:pipeline:general_information:signal_folder

Last update: **2026/09/09 20:57**

