

# Signal\_Names

/Signal\_Names or /SIGNAL\_NAMES The names of the signals used by the Pipeline Command.

The + symbol is used as a delimiter for multiple entries.

For example, to specify the three signals RFT1, RFT2, and RFT3

/Signal\_Names=RFT1+RFT2+RFT3 Many commands allow Signal\_Names to be blank

/Signal\_Names= In the cases that this is allowed, this indicates that Visual3D should use all signals of the specified [Signal\\_Types](#)

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:general\\_information:signal\\_names&rev=1720793066](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:general_information:signal_names&rev=1720793066)

Last update: 2024/07/12 14:04

