

An exception to this rule are [joint angles](#), which are represented in degrees.

Units_of_Measurement

Distance= m Velocity= m/s Acceleration= m/s² Mass= Kg Force= N Energy = J Power = W Moment
(Torque) = Nm Joint Angle = degrees Joint Angular Velocity = degrees/s Volume = m³

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:model_based_data_commands:units_of_measurement&rev=1721158234

Last update: 2024/07/16 19:30

