

calibration\_file or calibration\_file this command parameter refers to the link model.

if the parameter is empty, the *current* or *active* link model is used for the command.

the parameter should include the full pathname for the model or it may include wildcards.

note that the link model has the same name as the standing trial with which it is associated.

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:model\\_commands:calibration\\_file&rev=1718801522](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:model_commands:calibration_file&rev=1718801522)

Last update: 2024/06/19 12:52

