

# Set Segment Draw Options

## Overview

This command allows the user to...

## Pipeline Command

The command can be found in the **Pipeline Workshop** under **Model Building** as so:


```
Set_Segment_Draw_Options
/CALIBRATION_FILE=
/SEGMENT_NAME=
! /REMOVE_CURRENT_SHAPES=TRUE
! /DRAW_PRIMITIVE_SHAPE=TRUE
! /PRIMITIVE_SHAPES=FILLED_CIRCLE
! /PRIMITIVE_COLORS=BLUE
! /USE_TRANSPARENCY=TRUE
! /PRIMITIVE_LOCATIONS=VECTOR(0,0,0)
! /PRIMITIVE_SCALES=VECTOR(0.4,0.4,0)
! /PRIMITIVE_ROTATIONS=VECTOR(0, 0, 0)
! /PRIMITIVE_EVENTS=
;
```

## Command Parameters

|                         |                                                     |
|-------------------------|-----------------------------------------------------|
| <b>CALIBRATION_FILE</b> | The name of the hybrid model. Use the full pathname |
| <b>SEGMENT_NAME</b>     | The name of the desired segment                     |

## Dialog

The command can be edited in a text editor or in a dialog form. To edit in the dialog pop up form either click on the Edit button in the pipeline workshop or double-click on the pipeline command. The



dialog is shown below.

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:model\\_commands:set\\_segment\\_draw\\_options](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:model_commands:set_segment_draw_options)

Last update: 2026/09/21 18:17

