

Add Password Protection To Folder

The `Add_Password_Protection_To_Folder` command allows users to encrypt folders in the [data tree](#) with password protection.

NOTE 1: We strongly recommend that users DO NOT store data in their workspace that requires encryption, but some organizational policies require certain [TEXT_DATA](#) fields to be encrypted.

NOTE 2: There is no mechanism for recovering forgotten passwords. If the password used to encrypt a folder is lost then that data will no longer be readable.

Pipeline Command

The `Add_Password_Protection_To_Folder` command can be found in the pipeline workshop's Folder Encryption folder.

```
Add_Password_Protection_To_Folder
/SIGNAL_TYPES=
/SIGNAL_FOLDER=ORIGINAL
/Password=
;
```

Command Parameters

Parameter	Description
Signal Types	The type of signal whose folder will be encrypted.
Signal Folder	The signal folder to be encrypted.
Password	The password to be used. NOTE: Lost passwords cannot be recovered.

Notes

NOTE 1: We strongly recommend that users DO NOT store data in their workspace that requires encryption.

NOTE 2: There is no mechanism for recovering forgotten passwords. If the password used to encrypt a folder is lost then that data will no longer be readable.

See Also

Related pipeline commands include:

- [Remove_Password_Protection_From_Folder](#)
- [Enter_Password_To_View_Folder](#)

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:add_password_protection_to_folder

Last update: **2026/09/10 17:59**

