

Encrypt a folder with password protection.

Often used with [TEXT_DATA](#)

Add_Password_Protection_To_Folder /SIGNAL_TYPES= /SIGNAL_FOLDER=ORIGINAL /Password= ;

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:add_password_protection_to_folder&rev=1718648510

Last update: **2024/06/17 18:21**

