

Add_Password_Protection_To_Folder

Encrypt a folder with password protection.

Often used with [TEXT_DATA](#)

Add_Password_Protection_To_Folder /SIGNAL_TYPES= /SIGNAL_FOLDER=ORIGINAL /Password= ;

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:add_password_protection_to_folder&rev=1720027771

Last update: **2024/07/03 17:29**

