

Add Password Protection To Folder

Encrypt a folder with password protection.

Often used with [TEXT_DATA](#)

Add_Password_Protection_To_Folder /SIGNAL_TYPES= /SIGNAL_FOLDER=ORIGINAL /Password= ;

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:add_password_protection_to_folder&rev=1721229802

Last update: **2024/07/17 15:23**

