

this_is_a_real_time_command_that_can_be_used_to_draw_an_arbitrary_line_in_a_workspace. there_is_an_option_to_get_this_working_in_the_post_processing_signals_and_events_tab. example_signals_and_events_biofeedbackprocessingpipeline.png - open_a_cmz - signals_and_events_tab - go_to_view - view_biofeedback_post_processing - the_post_processing_dialog_should_appear_on_the_right_side_of_the_signals_and_events_tab - in_the_pipeline_dialog_under_rt_and_post_processing_biofeedback_select_animated_line_and_move_to_main_pipeline - click_edit_and_identify_which_signals_you_would_like_draw_a_line_between_in_the_command_code_animated_line_signal_types_target_target_signal_folder_original_original_signal_names_clav_cv7_line_width_2.0_line_color

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:animated_line&rev=1718801532

Last update: 2024/06/19 12:52

