

text_data_08_09_2017_by_wss_speed_0.4_a_forward_slash_is_a_special_character_for_pipeline_commands_that_indicates_that_a_parameter_name_follows_immediately_afterwards._if_an_equals_sign_follows_visual3d_confirms_that_a_parameter_was_specified._two_equals_signs_however_does_confirm_that_a_parameter_was_specified._visual3d_searches_for_an_equals_sign_after_each_slash_because_this_logically_indicates_that_the_user_has_specified_a_parameter._text_data_is_assigned_the_string_08_09_visual3d_finds_an_equals_sign_after_the_next_slash_so_interprets_the_rest_of_the_string_as_a_new_parameter_2017_by_mjt_using_modify_speed_metric_ver01_speed_0.4_it_doesn't_do_anything_with_this_parameter_because_it_doesn't_recognize_it._solution_remove_the_sign_text_data_08_09_2017_by_wss_speed_0.4_or_specify_two_equals_signs_text_data_08_09_2017_by_wss_speed_0.4}}}}}}}

From:
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:create_text_data&rev=1718801533

Last update: **2024/06/19 12:52**

