

Remove Password Protection From Folder

The Remove_Password_Protection_To_Folder command allows users to remove password encryption from [data tree folders](#).

Pipeline Command

The Remove_Password_Protection_To_Folder command can be found in the pipeline workshop's Folder Encryption folder.

```
Remove_Password_Protection_From_Folder
/SIGNAL_TYPES=
/SIGNAL_FOLDER=ORIGINAL
/PASSWORD=
;
```

Command Parameters

Parameter	Description
Signal Types	The type of signal whose folder will have password protection removed.
Signal Folder	The signal folder to have password protection removed.
Password	The encryption password.

See Also

Related pipeline commands include:

- [Add_Password_Protection_From_Folder](#)
- [Enter_Password_To_View_Folder](#)

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:remove_password_protection_from_folder&rev=1789067294

Last update:

2026/09/10 19:08

