

Set Animation Draw Size

Deprecated command in Version 6.

See: [Set_Animation_Draw_Properties](#)

Set_Animation_Draw_Size /SIGNAL_TYPES= ! /SIGNAL_NAMES= ! /SIGNAL_FOLDER=ORIGINAL !
/DRAW_SIZE=MEDIUM ;

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:set_animation_draw_size

Last update: **2024/07/17 15:46**

