

Set_Use_Processed_Targets

Sets a flag that causes Visual3D to use the signals in the TARGET PROCESSED folder for computing Kinetics and Kinematics.

Pipeline Command

Set_Use_Processed_Targets ! /USE_PROCESSED=TRUE ; see also [Using Processed Signals](#)

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:other_commands:set_use_processed_targets&rev=1720791051

Last update: 2024/07/12 13:30

