

Get a count of the number of parameter values in a parameter.

**Get\_Pipeline\_Parameter\_Count** /PARAMETER\_NAME= /PARAMETER\_VALUE= ;

## Example1

Create a parameter with three values, then retrieve the count

**Set\_Pipeline\_Parameter** /PARAMETER\_NAME=TEST\_PARAMETER  
/PARAMETER\_VALUE=ONE+TWO+THREE ; **Get\_Pipeline\_Parameter\_Count**  
/PARAMETER\_NAME=COUNT /PARAMETER\_VALUE=::TEST\_PARAMETER ; The results from executing these two command is:

---

**Set\_Pipeline\_Parameter** /PARAMETER\_NAME=TEST\_PARAMETER  
/PARAMETER\_VALUE=ONE+TWO+THREE ! /PARAMETER\_VALUE\_SEARCH\_FOR= !  
/PARAMETER\_VALUE\_REPLACE\_WITH= ! /PARAMETER\_VALUE\_PREFIX= !  
/PARAMETER\_VALUE\_APPEND= ; Pipeline Parameter : ::TEST\_PARAMETER = ONE+TWO+THREE

---

**Get\_Pipeline\_Parameter\_Count** /PARAMETER\_NAME=COUNT  
/PARAMETER\_VALUE=ONE+TWO+THREE ; ::COUNT = 3

## Example2

**Set\_Pipeline\_Parameter\_To\_List\_Of\_Signal\_Names** /PARAMETER\_NAME=TARGETS  
/SIGNAL\_TYPE=TARGET SIGNAL\_FOLDER=ORIGINAL ; **Get\_Pipeline\_Parameter\_Count**  
/PARAMETER\_NAME=COUNT /PARAMETER\_VALUE=::TARGETS ; For this specific test file the results are:

---

**Set\_Pipeline\_Parameter\_To\_List\_Of\_Signal\_Names** /PARAMETER\_NAME=TARGETS  
/SIGNAL\_TYPE=TARGET ! /SIGNAL\_FOLDER=ORIGINAL ; TARGETS= LAS LPP LPS LTH1 LTH2 LTH3 LTH4  
RAS RPP RPS RTH1 RTH2 RTH3 RTH4

---

**Get\_Pipeline\_Parameter\_Count** /PARAMETER\_NAME=COUNT  
/PARAMETER\_VALUE=LAS+LPP+LPS+LTH1+LTH2+LTH3+LTH4+RAS+RPP+RPS+RTH1+RTH2+RTH3+  
RTH4 ; ::COUNT = 14

From:  
<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:pipeline\\_commands:get\\_pipeline\\_parameter\\_count&rev=1718648529](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:pipeline_commands:get_pipeline_parameter_count&rev=1718648529)

Last update: **2024/06/17 18:22**

