

description: allow accessing pipeline parameter “array values” by “index”.

```
get_pipeline_parameter_index_value
/parameter_name=
/parameter_value=
/index=
;
```

example

```
set_pipeline_parameter
/parameter_name=running_speed
/parameter_value=10 km/h+12 km/h+14 km/h
;

set_pipeline_parameter
/parameter_name=filenames
/parameter_value=filea + fileb + filec
;

for_each
/iteration_parameter_name=file
/iteration_parameter_count_name=increment
/items=filenames
;
    file_open
    /file_name=:file
    ;
    get_pipeline_parameter_index_value
    /parameter_name= scott
    /parameter_value=:running_speed
    /index=:increment
    ;

    assign_tags_to_files
    /motion_file_names=:file
    /tags=:scott
    ;
end_for_each
/iteration_parameter_name=file
;
```

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:pipeline_commands:get_pipeline_parameter_index_value&rev=1718801548

Last update: 2024/06/19 12:52

