

Put_To_Matlab

Overview

Uses DCOM to put data from one active file in the Visual3D Workspace into Matlab (e.g. with Matlab running concurrently).

Pipeline Command

```
Put_To_Matlab  
/Matlab_Item_Name= the name of the signal in Matlab  
/Signal_Types= the Visual3D Type  
/Signal_Names= the Visual3D Name  
/Signal_Folder= the Visual3D Folder  
;
```

In Matlab the signal will appear with the Matlab_Item_Name. This is necessary because Matlab does not allow multiple signals with the same name but with a different type.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:pipeline_commands:put_to_matlab

Last update: **2025/09/12 19:40**

