

Append a saved report layout template.

```
Append_Report_Template  
! /REPORT_TEMPLATE=  
;
```

When the command executes a new Pipeline Parameter is created automatically.

EXISTING_NUM_REPORT_PAGES

This parameter contains the number of pages in the previously existing report.

This allows users to recognize the new page numbers in the appended report template.

Retrieved from ""

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:append_report_template&rev=1718644091

Last update: 2024/06/17 17:08

