

append a saved report layout template.

```
append_report_template  
! /report_template=  
;
```

when the command executes a new pipeline parameter is created automatically.

existing_num_report_pages

this parameter contains the number of pages in the previously existing report.

this allows users to recognize the new page numbers in the appended report template.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:append_report_template&rev=1718801561

Last update: 2024/06/19 12:52

