

Append_Report_Template

Append a saved report layout template.

```
Append_Report_Template  
! /REPORT_TEMPLATE=  
;
```

When the command executes a new Pipeline Parameter is created automatically.

EXISTING_NUM_REPORT_PAGES

This parameter contains the number of pages in the previously existing report.

This allows users to recognize the new page numbers in the appended report template.

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:append_report_template&rev=1721158283

Last update: **2024/07/16 19:31**

