

Export a metric table from the report to a text file.

```
Export_Metric_Table
/PAGE_NUMBER=
/COLUMN_NUMBER=
/ROW_NUMBER=
/EXPORT_FILENAME=
;
```

Retrieved from ""

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_metric_table&rev=1718386175

Last update: 2024/06/14 17:29

