

Export a metric table from the report to a text file.

```
Export_Metric_Table  
/PAGE_NUMBER=  
/COLUMN_NUMBER=  
/ROW_NUMBER=  
/EXPORT_FILENAME=  
;
```

From:

<https://has-motion.com/wiki/> - **Wiki Documentation Site**

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_metric_table&rev=1718648548

Last update: **2024/06/17 18:22**

