

Export_Metric_Table

Export a metric table from the report to a text file.

```
Export_Metric_Table
/PAGE_NUMBER=
/COLUMN_NUMBER=
/ROW_NUMBER=
/EXPORT_FILENAME=
;
```

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_metric_table&rev=1719433722

Last update:

2024/06/26 20:28

