

# Export\_Metric\_Table

Export a metric table from the report to a text file.

```
Export_Metric_Table
/PAGE_NUMBER=
/COLUMN_NUMBER=
/ROW_NUMBER=
/EXPORT_FILENAME=
;
```

From:  
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:  
[https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report\\_commands:export\\_metric\\_table&rev=1720027790](https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_metric_table&rev=1720027790)

Last update: 2024/07/03 17:29

