

Export_Metric_Table

Export a metric table from the report to a text file.

```
Export_Metric_Table
/PAGE_NUMBER=
/COLUMN_NUMBER=
/ROW_NUMBER=
/EXPORT_FILENAME=
;
```

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_metric_table&rev=1721149204

Last update: 2024/07/16 17:00

