

Merge a graph of signals into an ODT document.

[Export_Signals_To_ODT_Graph.png](#)

```
Export_Signals_To_ODT_Graph
! /FILE_NAME=
/GRAPH_NAME=
! /OWNERS=
/XSIGNAL_TYPE=
! /XSIGNAL_FOLDER=ORIGINAL
! /XSIGNAL_NAME=
! /XSIGNAL_COMPONENTS=
! /XSIGNAL_RANGE=NORMALIZE
! /XSIGNAL_MIN=
! /XSIGNAL_MAX=
/YSIGNAL_TYPE=
! /YSIGNAL_FOLDER=ORIGINAL
! /YSIGNAL_NAME=
! /YSIGNAL_COMPONENTS=
! /YSIGNAL_RANGE=MAXMIN
! /YSIGNAL_MIN=
! /YSIGNAL_MAX=
! /LINE_COLOR=#000000
! /LINE_STYLE=SOLID
! /LINE_LEGEND=
! /EVENT_SEQUENCE=
! /EXCLUDE_EVENTS=
! /NORMALIZE_DATA=FALSE
! /NORMALIZE_POINTS=101
! /EXPORT_MEAN_AND_STD_DEV=FALSE
;
```

From:
<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:
https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_signals_to_odt_graph&rev=1718648549

Last update: 2024/06/17 18:22

