

Export_Signals_To_ODT_Graph

Merge a graph of signals into an ODT document.



```
Export_Signals_To_ODT_Graph
! /FILE_NAME=
/GRAPH_NAME=
! /OWNERS=
/XSIGNAL_TYPE=
! /XSIGNAL_FOLDER=ORIGINAL
! /XSIGNAL_NAME=
! /XSIGNAL_COMPONENTS=
! /XSIGNAL_RANGE=NORMALIZE
! /XSIGNAL_MIN=
! /XSIGNAL_MAX=
/Y SIGNAL_TYPE=
! /Y SIGNAL_FOLDER=ORIGINAL
! /Y SIGNAL_NAME=
! /Y SIGNAL_COMPONENTS=
! /Y SIGNAL_RANGE=MAXMIN
! /Y SIGNAL_MIN=
! /Y SIGNAL_MAX=
! /LINE_COLOR=#000000
! /LINE_STYLE=SOLID
! /LINE_LEGEND=
! /EVENT_SEQUENCE=
! /EXCLUDE_EVENTS=
! /NORMALIZE_DATA=FALSE
! /NORMALIZE_POINTS=101
! /EXPORT_MEAN_AND_STD_DEV=FALSE
;
```

From:

<https://has-motion.com/wiki/> - Wiki Documentation Site

Permanent link:

https://has-motion.com/wiki/doku.php?id=visual3d:documentation:pipeline:report_commands:export_signals_to_odt_graph&rev=1721158286

Last update: 2024/07/16 19:31

